

Electron Spin Resonance Spectrometer

KS/PHY/ELECPRESSPEC

Exp-Determination of Lande-g-factor

Required Components For the experimental Set-Up

Complete Experimental Set up

Experimental Kit

With Helmholtz coils (152mm dia with 500 turns) 180x180x185mm separate coil arrangement

Pickup coil with 40 turns and DPPH sample

Built-in Digital frequency meter

Complete unit housed in metal cabinet of 385x210x165mm with effective electromagnetic shielding

Proper cables RG-59 for connection with oscilloscope included with the setup

➤ FET based marginal R.F. Oscillator

➤ Digital display of frequency

➤ Excellent peaks display

➤ Digital display of Helmholtz Coil Current

➤ Compatible with general purpose CRO in X-Y mode

➤ Advantages and Limitations of our Spectrometer:

➤ The instrument is basically designed for postgraduate laboratories keeping in view their requirements and limitations.

➤ The observation of ESR at low magnetic fields and consequently in radio-frequency region makes its instrumentation and working a lot simple and within the reach of a postgraduate student. Good resonance peaks can be obtained as a class room exercise.

The spectrometer is complete in all respects including a sample DPPH.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Digital Gauss Meter

KS/PHY/DIGGAUSSMTR

Purpose

A Digital Gauss Meter measures magnetic flux density (magnetic field strength) in gauss (G) or tesla (T). It is used in physics labs, material testing, magnet characterization, and electromagnetic device analysis.

Principle

Most digital Gauss meters operate using a Hall-effect sensor. When the sensor is placed in a magnetic field, a voltage proportional to the field intensity is generated. The meter amplifies, conditions, and digitizes this voltage to display magnetic field strength.

Key Components

Hall probe

- Axial probe for fields parallel to the probe axis

- Transverse probe for fields perpendicular to the probe axis

Probe handle and cable

- Shielded to minimize noise

Control and display unit

- High-gain amplifier

- Analog-to-digital converter

- LCD/LED digital readout

- Range selector

Calibration magnet or zero-field chamber

- For zeroing and calibration



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Digital Gauss Meter

KS/PHY/DIGGAUSSMTR

Power supply

Internal battery or DC adapter

List of Experiment

Measurement of magnetic field

Digital Gauss meter with GaAs Transverse Probe

Range: 0-2KG, 0-20KG and

7 segment LED display with accuracy 0.5%

Complete unit housed in metal cabinet of 260x200x110mm with effective electromagnetic shielding

Cabinet shall be powder coated with plastic moulding on edges

Range : 0-20 K.Gauss/ 0-20 K Gauss

Accuracy 1% of reading

Display : LCD display (16X2) with backlight

Microcontroller based



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Constant Current Power Supply

KS/PHY/CONSCURRNTPS

Purpose

A constant current power supply delivers a fixed, regulated current to a load regardless of changes in load resistance or input voltage. It is used for LED driving, battery charging, sensor biasing, and laboratory experiments requiring stable current.

Principle

Current regulation is achieved by sensing the load current through a resistor and controlling the output via a feedback loop so that:

$$I_{out} = V_{ref} / R_{sense}$$

Any change in load voltage or resistance is compensated by adjusting the output voltage to keep the current constant.

Core Components

Current-sense resistor

Low-value resistor for accurate current measurement

Error amplifier or op-amp

Compares sensed voltage with reference

Reference voltage source

Provides a stable set-point

Series pass element

Power transistor or MOSFET

Adjusts voltage to maintain constant current



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Constant Current Power Supply

KS/PHY/CONSCURRNTPS

DC supply input

Often 12 V, 24 V, or adjustable

Protection circuitry

Thermal, over-voltage, and short-circuit protection

Output terminals with current indicator

Digital meter or LED display

List of Experiments

To provide constant current to an electromagnet

Constant Current Power Supply (transistor based) for an electromagnet (0-4A current) continuously variable.

Complete unit housed in metal cabinet of 390x280x160 mm with effective electromagnetic shielding

Cabinet shall be powder coated with plastic moulding on edges.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Photoconductivity

KS/PHY/PHOTOCONDCVTY

List of Experiments

To study the current vs voltage characteristics of CdS photoresistors at constant irradiance

To measure the photo current as a function of irradiance at constant voltage

OPTICAL BENCH

Material: Aluminium alloy

Type: Hexagonal Section

Scale: 0 - 100 cm

Least count: 1 mm

POWER SUPPLY 0-16V

Voltage: 0-16 Volt DC continuously variable & stabilized

Ripple: Less than 25mV

Overload protection: Current limiting

Current: 1 Amp. continuously variable from 10% to full rating

Display: Two separate displays (3 digit LED) are provided to monitor the output voltage and load current continuously

Working voltage: 230V AC, 50 Hz single phase

PHOTORESISTOR LDR IN MOUNT

Photo resistor: CdS (Cadmium sulphide)

Aperture: 10 mm, Clear

Mounting rod: 10 mm dia.

Connection: 4mm safety terminals

Working voltage: 0-16V DC



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Photoconductivity

KS/PHY/PHOTOCONDCVITY

CONVEX LENS IN HOLDER

Focal Length: 100mm

Diameter of Lens: 50 mm

Frame Diameter: 130 mm to avoids scattering of lights

Rod Diameter: 10 mm

POWER SUPPLY 12V AC/DC

Output: 2, 3, 4, 5, 6, 8, 10 & 12 VAC full wave rectified, unsmoothed & unregulated DC

Overload: Resettable thermal trip.

Input: 230 V AC, 50 Hz

LAMP HOUSING

Lamp: 12V, 21V

Lens: Spherical condenser, to and fro adjustment

Connection: 4mm plug lead

Mounting rod: 10 mm dia.

Housing Aluminium, Heat ventilation arrangement



Photoconductivity Experimental Set up: Setup on a aluminium profile type optical bench with rides, halogen light source variable power with housing. Mounted photo conductor with BNC and Cable, DC regulated power supply with meter, lens with lens mount. Filters in filter when in visible wavelength (400-700nm) with FWHM 10nm and various accessories like DMM and cables.

Learning Objectives/experiments

- < To study the current vs voltage characteristics of CdS photo-resistor at constant irradiance.
- < To measure the photo-current as a function of the irradiance at constant voltage.

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Faraday Effect with Diode Laser

KS/PHY/FARADYEFF

List of Experiments

Exp-1: Observing the rotation of the polarization plane when polarized monochromatic light passes through flint glass under influence of a magnetic field.

Exp-2: Determining Verdet's constant from the relation between rotation angle and magnetic flux.

Exp-3: Verification of the relationship between Verdet's constant and wavelength

Flint Glass: DEDF glass, size 15x15x15mm

HC Power supply 0-30V DC, 10A, 3 1/2 Display LED

Digital Gauss meter

Range 2KG & 20 KG,

Resolution 1G & 10G,

Sensor: In As


Electromagnet coil: 500 T coil, 7Amp (Max), 185WG, Cu,

Electromagnet Core: Ferromagnetic, 40x40mm cross section,

Optical bench

Al exclusion,

L= 1mtr,



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Faraday Effect with Diode Laser

KS/PHY/FARADYEFF

Triangular shape with end stop, 4 SS leveling screw

Transversal slider: Al, 25-0-25mm, LC 0.01mm

Fixed slider: Al, Base width 50mm, Height 115mm

Diode laser

Wavelength 532nm, 5mW, Class II

Scope of delivery

Flint glass with Holder

Electromagnet unit with slider

Polarizer/Analyzer

Object screen

Optical bench, 1 m

Fixed slider

Power supply (0-30V DC, 10A)

Digital gauss meter with probe

Hall probe holder

Convex lens ($F=10\text{cm}$)

Cylindrical base

Green Diode laser

Power Adapter



OPTICAL BENCH TRIANGULAR Material : Aluminum extrusion

Type : Triangular shape Scale : 0-100cm

Least count : 1mm

This optical bench is rigid, heavy, stable and long lasting.

It has

four levelling screw and flexible feets. **He-Ne LASER** Wavelength :

632.8 nm Working current : 4mA ~ 6mA

Output power : > 2mW Working time : > 8 hrs.

Working voltage : AC 220 V $\pm$ 22 V Input Power : <2 W

Dimension : 300x62x82 mm Weight : 1.5 kg (approx.)

DIODE LASER

Peak wavelength : 635nm Operating voltage : 5V DC

Operating current : 250mA Optical power : 0.4-0.8mW

Laser product : Class II Operating temp. : 0 - 40°C

Storage temp. : -10 to 50°C

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Study of pn Junction & Energy Band Gap

KS/PHY/PNJN&ENGYBNGP

List of Experiment

Exp-1 Determination of reverse saturation current.

Exp-2 Study of Energy Band Gap of p-n Junction

Exp-3 Study of Junction capacitance

P-N JUNCTION SETUP

Selector Switch: V-I and V-T experiment, Bias & Junction

Selector Switch at V-I position/Junction

Voltmeter Display: 3 1/2 digit, 7 segment LED, auto polarity

Voltage Range: 0.000-1.999V

Current Display: 3 1/2 digit, 7 segment LED

Current Range: 0-20mA

Selector Switch at V-T position/Junction

Voltage Display 3 1/2 digit, 7 segment LED

Voltage Range: 0.000-1.999V

Temperature Display: 3 1/2 digit, 7 segment LED

Temperature Range: 273 K to 353 K

CRO in Bias Position

Frequency: 5 KHz & 20 KHz

Voltage: 220 mV (p-p)


Output Connector: 3 Pin, DIN type

Voltage Range: 0.00-10.00V

Oven Connector: 5 Pin, DIN type

Diode & Transistor: 4 mm safety socket

Input Voltage: 220 V, 50 Hz AC



OVEN WITH TEMPERATURE SENSOR

Heating Element: 35 ohm

Oven Connector 5 Pin, DIN type

Ambient Temperature 353

Temperature Sensor: Pt 100

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Kerr Effect Experiment Kit

KS/PHY/KERREFFEXP

Purpose

A Kerr Effect Experiment Kit is used to study the **Electro-Optic Kerr Effect**, where an applied electric field induces birefringence in an optical medium. This results in rotation of light polarization, allowing measurement of Kerr constants and modulation of light intensity.

Principle

When a strong electric field E is applied across a Kerr medium (commonly nitrobenzene), the refractive index changes proportionally to E^2 :

$$\Delta n = K E^2$$

This induced birefringence converts linearly polarized light into elliptically polarized light. Using crossed polarizers, the intensity change is measured as a function of voltage.

Exp-1 To demonstrate the Kerr effect in nitrobenzene solution (only observation)

Exp-2 To measure the light intensity as a function of voltage across the Kerr cell using photodetector

Kerr Cell

Electrode gap: 1mm

Voltage limit: 5 kV, DC (max)

Dimensions: 50x50x20 mm

Polarizer/Analyzer


Angle: Adjustable (0 - 90 deg)

Aperture: 21mm dia

Frame 130mm dia Blackrood, to avoid scattering of light

Rod: 10 mm dia

Least count: 1 deg



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Kerr Effect Experiment Kit

KS/PHY/KERREFFEXP

Halogen Light Source

Halogen bulb :12V, 50W

Operating voltage: 12V, 5A

Safety sockets: 4 mm

Rod dia: 10mm

Object screen


High voltage power supply

Input voltage: 230V AC 50Hz

Output Voltage: 0- 6KV DC: 6.3V AC/2 Amp

Display: 3 1/2 digit LED

Current limit: 50 μ A & 2 mA(Max.)



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Kerr Effect Experiment Kit

KS/PHY/KERREFFEXP

Power Supply

Output: 2,3,4,5,6,8,10 & 12VAC full wave rectified, unsmoothed & unregulated D.C.

Overload: Resettable thermal trip.

Input: 230 VAC, 50 Hz

Prism Table

Optical bench triangular

Material: Aluminum extrusion

Type: Triangular shape

Scale: 0-100cm

Least count: 1mm

This optical bench is rigid, heavy, stable and long lasting. It has four levelling screw and flexible feets

Convex lens in holder

Set of six filter

Flexible plug lead

Transversal slider

Pin hole photodetector

Digital Multimeter

Power cord

Fixed slider

Fixed optical slider



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Zeeman Effect Experiment Kit

KS/PHY/ZEEMAN EFF

Purpose

The Zeeman Effect Experiment Kit is used to study the splitting of atomic spectral lines when an external magnetic field is applied. This experiment verifies electron orbital magnetic moment theory and allows determination of the Landé g-factor.

Principle

When an atom is placed in a magnetic field B , its energy levels split according to:

$\Delta E = m_l \mu_B B$, $\Delta E = m_l \mu_B B$ This causes a single spectral line to split into multiple components (normal or anomalous Zeeman effect). Using an interferometric or high-resolution spectroscopic setup, these split components are observed and measured.

Experiment

EXP 1. To observe the Zeeman splitting of the green (546.1nm) mercury line using Fabry-Perot etalon for normal

EXP 2. transverse and longitudinal configuration


EXP 3. Determination the polarization as of the triplet components in transverse configuration

EXP 4. Determining the polarization state of the doublet components in longitudinal configuration

Electromagnets

Ferromagnetic core, Coils 500 T. Current 5A (max), 4 mm safety socket

U Core 150 x 130 mm (LxH), 40 x 40 mm, cross section, I Core L=150mm



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Zeeman Effect Experiment Kit

KS/PHY/ZEEMANEFF

Fabry Perot Etalon

Mirror optics $\lambda/20$, Mirror gap 3mm, Filter 532nm, Green Interference Filter, Clear View 40mm dia, Rod 10 mm dia

CMOS Camera

1.3 MP, USB Interface

Power Supply


0-40V DC, 0-10A, voltage resolution 0.1V, display 3 1/2 digit LED, Current resolution 0.1A, primary fuse 8A

Optical bench

Al exclusion, 1-1.5 mm with Scale, Triangular shape with endpoint leveling (steel screw)

Optics holder

Al Dia 130 x 3 mm ($\Phi \times T$), H 230 mm



Kalpna Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Zeeman Effect Experiment Kit

KS/PHY/ZEEMANEFF

Traversal slider

25-0-23 LC 0.01mm

Scope of delivery

Optical Bench briangular

Transverse slider

Convex Lens in Lens in Holder, FL. 150 mm

Convex Lens in holder, FL 50mm

Polarizer fiber

Quarter wave filter

Electromagnet Coil

Mercury Lamp with Power Supply

Mercury Lamp holder

Fabry-Perot etalon

CMOS Camera with Software


Digital Gauss meter

Green Filter with holder

Micrometer Eyepiece

Fixed slider

Desktop Computer



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Purpose

To determine the frequency of a vibrating source or the tension in a string using standing waves formed on a stretched string driven by a vibrator.

Principle

A string under tension T and linear mass density μ supports standing waves.

The wavelength is related to the string length and mode number.

Wave velocity:

$$v = \sqrt{T/\mu}$$

Frequency:

$$f = v/\lambda$$

Modes for a string fixed at both ends:

$$\lambda_n = 2L/n$$

Experiment

EXP 1. To determine the frequency of an electrically maintained tuning fork by Melde's Experiment & verify $\lambda^2 \propto T$

Heavy steel fork

Heavy cast iron base

Electromagnet

Weight box

Voltage source 1.5 V – 12 V/ 3A

Pulley with clamp

Reel of thread

Meter scale of length 1m

Scale pan

Aim

To determine the frequency of an electrically maintained tuning fork by,

1. Transverse mode of vibration
2. Longitudinal mode of vibration

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com

Melde's Experiment

KS/PHY/MELDEEXP



Michelson Interferometer Kit

KS/PHY/MICHELSNINTRFE

Purpose

To observe interference fringes and perform precision measurements such as wavelength determination, refractive index variation, and small displacement measurement.

Principle

Light from a coherent source is split into two beams by a partially silvered **beam splitter**. These reflect off two mirrors (one movable) and recombine to form interference fringes. Moving one mirror by distance d shifts fringes: $2d = m\lambda$. To determine the wavelength of monochromatic light (using sodium light source)

EXP-2. To determine the wavelength of monochromatic light (using He - Ne laser).

Moveable mirror: Range 10-0-10mm

Course knob Least count 0.01mm

Micrometer fine: Least count 0.0001mm

Beam splitter: Rectangular, size 50 x 29 x 7mm. R/T% 50/50, Flatness $\lambda/8$ at (632 nm)

Compensator Rectangular, size 50 x 29 x 7mm, Flatness $\lambda/8$ at (632 nm)

Mirror Circular: 2 no, dimension 30x10mm

Dimension: 290 X 210 X 230mm (LxWxH)

Telescope : Tube 90 X 25mm (L x Φ) , X 10 eye piece


Sodium source: MS housing, 200 X 85mm (L x Φ) slit (L x W) 20 x 16 mm

Laser: He-Ne Laser (632nm), Power 2 mW

Beam Expander with holder

object screen with cylinder base

Lab jack



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Purpose

To verify the quantized nature of atomic energy levels by observing inelastic collisions between electrons and gas atoms (typically mercury or neon).

Frank Hertz Experiment

KS/PHY/FRANKHRTZ

Principle

Electrons emitted from a heated cathode are accelerated through a low-pressure gas. When electron kinetic energy reaches specific threshold values equal to atomic excitation energies, they undergo inelastic collisions, losing energy. This causes periodic drops in the collected anode current.

Key relation:

$$eV = \Delta E$$

The spacing between successive minima/maxima in the I-V curve equals the excitation energy of the gas.

Frank hertz experiment (Neon tube) with digital data logger & software

EXP 1. To observe the neon spectral bands formation in Frank Hertz tube

EXP 2. To record the Frank Hertz characteristics curve for neon

EXP 3. To study the effect of filament voltage and anode plate on the characteristics plate

Frank Hertz unit: mode manual or data logging

Voltage: 1.2-12 V, Display: 3 1/2 digit LED screen

Grid(G) voltage : 0 -80

Control Grid(G) voltage : 0 -10V

Filament Voltage: 0 -9 V

Neon Tube:

Dia: 130 x 26 mm

Filament Voltage: 6 -9 V

Control voltage: 9 V

Accelerating Voltage Max : 80 V

Counter Voltage: 1.2 - 10 V



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Frank Hertz Experiment

KS/PHY/FRANKHRTZ

Data logger

I/P channels 4, 6 PIN BT connector,

Resolution 12 bit ADC @100 ksps

Analog I/P 0-5V,

O/P channels 12 bit DAC

Analog O/P +/- 12 V max 10 mA, Current booster

1A, Power connector 3 PIN, DIN, USB cable

Power Supply: 220V, 50Hz AC output +/- 12 V, 3 PIN DIN

Voltage Sensor +/- 10V, 4mm socket, 6 PIN BT connector

Current Smuor: 100A, 4mm socket, 6 PIN BT connector

In Package

Franck-Hertz unit

Neon tube with mount

Power Cord

Data logger

Power supply for data logger

Current sensor +/- 100mA

Voltage sensor maor +/- 10V

Software CD disc, With Desktop Computer



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Automatic Spin Canting System & Accessories

KS/PHY/ASCS&ACCESS

Purpose

Used for precise measurement and control of spin canting angles in magnetic materials, nano-magnetic samples, thin films, and ferrites.

Provides automated rotation, magnetic field alignment, and data acquisition for magneto-optical and magnetization

studies. Programmable Speed Range: 50 -12,000 RPM

[based on a glass substrate of dimension 38 mm (L) x 25 mm (W) x thickness 1 mm (H)]

The and on a Glass Substrate Dimension 18 mm (L) X 25 mm (W) & Thickness 1 mm (H)

Programmable Acceleration: 40-6,000 RPM/s

[based on a glass substrate of dimension 38 mm (L) x 25 mm (W) x thickness 1 mm (H)]

Can be on a Glass Substrate of Dimension 18 mm (L) X 25 mm (W) & Thickness 1 mm (H)

< +/- 1% Error across the Full Speed Range

N₂ & other Inert Gas Purging Port

Programmable Controlling Duration: 1-9,999 sec./step

30% editable program

50 editable steps / program

Microcontroller controlled

PC connectivity through USB 2.0 (for PC mode)

warm up option

Calibration Option

Spill-drainage facility

High-speed DC Motor

Non Volatile Program Memory

user friendly firmware interface


Dimension: 291 mm (L) X 448 mm (W) X 243 mm (H)

Input & Controlling through Key-pad

PTFE coated Working Chamber of 8 " Diameter

Integrated Power On/Off Switch with Indicator

Integrated Lid Protection Switch with On/Off Indicator



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Automatic Spin Canting System & Accessories

KS/PHY/ASCS&ACCESS

Integrated Mechanical Lid Inter-lock with Lid-open Sensor

Integrated Vacuum Release Switch with On/Off Indicator

Transparent Photo-resist Safety Lid over the Working Chamber

In-situ Sample Dispensing Port on the Transparent Photo-resist Safety Lid

Real-time Display of RPM, Timing & Program Status on Graphical LCD Console

Real-time Display of R.P.M. during Warm-up & Calibration on Graphical LCD Console

Real-time Display of R.P.M. vs. Timing Graph during Program Run-time on Graphical LCD Console

Power Input: Universal

Advanced Oil-free Vacuum Pump

Power Input: Indian Standard

Power: 80 W (Maximum)

Overflow Protection available

Vacuum-735 mm/Hg. (Maximum)

Flow Rate: 6-10 l/min

Motor Rotation 1,450 RPM.

Capacitance 3 micro farad 450V

Horse Power : 1/6 HP

Net Weight: 5.2kg

Port Thread: 0.3125"


Noise Level: 52 dB

Current: 0.4 A (Maximum)

Moisture Trap available

Vacuum Regular available

Dimensions: 310 mm (L) X 135mm(W) X 204(H)



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Viscosity by Stoke's Method

KS/PHY/VISCOSTYSTOKE

Purpose

To determine the coefficient of viscosity of a viscous liquid (usually glycerin) by measuring the terminal velocity of a falling spherical ball.

Principle

A small sphere falling through a viscous liquid experiences:

Gravitational force

Buoyant force

Viscous drag force (Stokes' law)

To determine the coefficient of viscosity of glycerin by falling sphere method

To determine the density of given liquid using steel sphere, glass beads

To predict the fall time of different spheres of same material

Tube Stand

base (cast iron): 23 x 15 cm

Rod mild steel: 110 cm

Cylindrical Tube

Length 100 cm (approx)

Internal diameter: 3.5 cm

External diameter: 4.0 cm

Volume : 962 cm³

Measurement Unit

Mains Supply Voltage: 230V +/- 10%, 50Hz

Adaptor Output: 5V DC

Timer checking time: 5 sec

Time segments: 3

Steel Sphere

Diameter 0.2 cm to 0.5cm


Density: 7.85 gm/cm³

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Semiconductor Diode Laser with Power Supply

KS/PHY/SEMIDIODLASR

Purpose

A semiconductor diode laser with its dedicated power supply is used for optics, fiber-optic communication experiments, holography, interferometry, and opto-electronic characterization. It provides a stable, coherent, monochromatic light source.

Principle

A p-n junction diode with a direct bandgap material (GaAs, AlGaAs, InGaAs, etc.) emits coherent light via **stimulated emission** when forward-biased above threshold current. The diode is housed in a laser package with mirrors forming a Fabry-Perot cavity or distributed feedback (DFB) structure.

A compact, fully solid-state instrument produces an intense beam of light at a wavelength 635 - 670 nm for red color. This semiconductor diode laser is mounted in a simple mount. Supplied completely with power supply.

Operating voltage 3V

Optical power 3-5mW

Wavelength: 635nm

Optical Power 3-5mw

Operating Current <250mA

Operating Voltage 3V DC

Operating Temperature 15°C-35°C

Storage Temperature: 0°C-50°C

Spot Size Range: 300µm

Spot Size: < 6.5 mm

Lifetime >3000 hrs



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



IC Tester

KS/PHY/ICTESTR

Purpose

An IC Tester checks the functionality of digital and analog integrated circuits by applying predefined logic patterns or electrical conditions and comparing outputs with expected responses.

Types

Digital IC Tester

- Tests logic families: TTL, CMOS

- Supports ICs like gates, counters, flip-flops, decoders, multiplexers

Analog IC Tester

- Tests op-amps, comparators, voltage regulators

Universal IC Tester

- Programmable patterns

- Large IC library

- Automatic pin detection and fault diagnosis

40 Pin Digital IC Tester to test 54/74 series TTL, and 4000 & 4500 series CMOS ICs, tests more than 2000 kinds of devices

Self Diagnosis

Identifies Unknown model number of Devices

- Measures more than 2000 kinds of devices

test 54/74 series TTL

test 4000 & 4500 series CMOS

40 Pin capability

Control: 16 key test switch keypad with dual tone sound indication

Display: 6 Digit LED Display

Test Socket: 40 pin ZIP Socket

Library of ICs: TTL54 Series, TTL55 Series, TTL74 Series, TTL75 Series CMOS 14 Series, CMOS 40 Series, Optical Coupler Series, LED Display Series, RAM Series, SCM Series, CPU Peripheral Series

Auxiliary Power Supply: 220V AC 50Hz

Kalpna Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Purpose

To study the **linear electro-optic (Pockels) effect**, where an applied electric field changes the refractive index of a non-centrosymmetric crystal (e.g., KDP, LiNbO₃). Used to observe phase retardation, intensity modulation, and determine the electro-optic coefficient.

Principle

In a Pockels crystal:

$$\Delta n \propto E$$

Applying a voltage produces a phase shift between ordinary and extraordinary rays. With crossed polarizers, the light intensity varies with voltage. At half-wave voltage V_{π} , the induced phase shift is π .

To plot the graph and study the birefringence with respect to applied voltage in an electro optic crystal (Lithium niobate)

Polarizer Rotator with Mount

Resolution 2 deg

Rotation 360 degree

Polarizer with Precision Rotation

Material: Sheet Polarizer

Rotation: 360 deg


Fine Rotation : 15 degree

Resolution: 0.1 degree

Optical Rail

Length: 1000mm

Material: Black anodized Aluminum alloy



Pockel Effect Apparatus

KS/PHY/POCKELEFF

Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Pockel Effect Apparatus

KS/PHY/POCKELEFF

Kinetic Laser Mount

Material: Black anodized Aluminum alloy

Adjustment: Using 80 tpi lead screws

Adjustment Range: ± 3 degrees

Detector Mount with X-Translation

Material: Copper Black Anodized Aluminum Alloy

Resolution: 0.01mm


Travel: Micro meter control

Pockel Cell with Mount

Electrode Material: copper

Cavity width : 0 - 6 mm Lead Screw control

Cavity Size: 30x5x5 mm



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Pockel Effect Apparatus

KS/PHY/POCKELEFF

Electro-Optic Crystal

Material: Lithium niobate (Li NbO_3)

Dimension : 25x4x4mm

Switched Mode Power Supply

SMPS Input: 230v 50 Hz

Voltage: 0-2 KV

Diode Laser with Power supply (Red)

Wave length: 650 nm


Optical power: 3mW

Detector Output Measurement Unit

Sensor Type: Photo Transistor

Display: 7 segment, 3 1/2 digit

Range: 0-199 mill/micro



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Lees Disc Apparatus

KS/PHY/LEEDISCAPP

Purpose

Used to determine the thermal conductivity of a poor conductor (insulating material such as cardboard, felt, ebonite, or rubber).

Principle

At steady state, the heat conducted through the test specimen equals the heat lost from the upper disc by radiation and convection. Using temperature differences and heat flow rate, thermal conductivity KKK is calculated.

Main Components

EXP- To determine thermal conductivity of bad conductors (Glass, Nylon, Plywood, etc.) in the form of a disc using Lee's method

Required Components For the experimental Set-Up

Complete Experimental Set up

Disc Sample type: Cardboard, Glass, Plywood

Temperature Sensors: 0.5°C (Thermometer)

Disc Diameter: 111 mm

Disc Thickness: 2.8 mm

Components: Electric Oven for heating, heat chamber, Sand, rubber tube, water boiling container ect.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Searls's Method Apparatus

KS/PHY/SEARLMAPP

Purpose

To determine the **Young's modulus** of a wire by measuring its extension under load in a controlled setup.

Principle

A wire under tensile load extends proportionally to the applied force:

$$Y = FL / A\Delta L$$

Searle's method compares the extension of the test wire with a reference wire to eliminate errors due to frame movement and temperature effects.

To determine conductivity of good conductor

Required Components For the experimental Set-Up

Complete Experimental Set up

Copper Turing, Steam Chamber, Digital Scale, Thermometer, Beaker, Tubing, Maker, insulating Material, Vernier Caliper, Metal Rod.

Consists of a cylindrical copper rod about 300 mm long & 25 mm diameter with its one end surrounded by a steam jacket for heating it from a boiler & the other end is kept cool by a system of water flowing through a spiral tube soldered to that end. Through two sockets, two thermometers are inserted in the rods for finding the temperature at two points. Two more thermometers are inserted in the other end which is kept cool by the stream of water. Through this apparatus, the difference of temperature between the inflowing & outflowing water is determined. Fitted in teak wood case, packed with felt with removable front; with thermometers.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Optical Fiber Characterization Apparatus

KS/PHY/OPTFIBCHAR

Purpose

Used to measure key parameters of optical fibers such as numerical aperture, attenuation, bending loss, dispersion, and coupling efficiency. Exp. Numerical aperture measurement of multi mode fibre, measurement of bending loss in multi-mode fiber, Numerical aperture measurement of single mode fibre, Relative measurement of splice loss in multi mode fibre, Calculation of normalized Frequency or V- number of single mode fibre, Calculation of mode field diameter of single mode fibre.

Required Components For the experimental Set-Up

Complete Experimental Set up

Optical Rail

Length: 1000 mm

Material: Aluminum alloy

Finish: Black anodized

Quantity: 1 no.

XYZ Translation stage with mount:

Material: Black anodized aluminium alloy

Travel: Micrometer controlled

Resolution: 0.01mm

Quantity: 1 no

Fiber Check Holder:

Material: Aluminum alloy

Finish: Black anodized

Quantity: 1 no.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Optical Fiber Characterization Apparatus

KS/PHY/OPTFIBCHAR

Single Mode Optical Fiber

Material: Glass Fibre

Numerical Aperture: 0.11

Core Dia: 9 micron

Length: 3 meters

Quantity: 1 no=

Detector Output Measurement Unit

Sensor Type: Photo Resistor

Display 7 segments, 3 1/2 digit

Range 0-199 milli/micro ampere

Quantity: 1 no

Kinetic Laser Mount:

Material: Black anodized aluminium alloy

Adjustment: Using 80 tpi lead screws

Adjustment Range: +/- 4 degrees

Quantity: 1 no

Laser Fibre Coupler with Mount:

Magnification of object: 10 x

Positioning: lead screw controlled

Quantity: 1 no

Bending Loss Apparatus

Chuck material: Nylon

Step diameter 35, 45, 55, 65 mm

Quantity: 1 no



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Optical Fiber Characterization Apparatus

KS/PHY/OPTFIBCHAR

Multi Mode Optical Fibre

Material: Plastic fibre

Numerical Aperture: 0.5

Core Diameter: 735, 240 microns

Length: 5 meters

Quantity: 1 no each

Detector Mount

Material: Black anodized Aluminum alloy

Diameter: 30 mm

Quantity: 1 no

Fibre Holder with Angular Tilt

Material: Aluminum alloy

Finish: Black anodised

Resolution: 2 degree

Quantity: 1 no

Fibre Chuck:

Material: Black anodized Aluminum alloy

Diameter : 30 mm

Quantity: 1 no

Diode Laser with Power supply (Red)

Wave length: 650nm

Output power: 5mw

Quantity: 1no



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Apparatus for Opto Electronics Characterization

KS/PHY/OPTOELECHAR

Purpose

To characterize the electrical and optical properties of opto-electronic devices such as LEDs, laser diodes, photodiodes, phototransistors, LDRs, and solar cells.

Light Sources

LED sources (various wavelengths)

Laser diode modules

Halogen or tungsten lamps

Monochromator (optional for wavelength-selective illumination)

2. Photodetectors

Photodiodes (PIN, avalanche)

Phototransistors

LDRs

Solar cells

Photoconductive cells (CdS, CdSe)

3. Optical Components

Collimating lenses

Optical filters

Diffusers

Beam splitters

Optical fiber patch cords

Adjustable apertures

4. Mechanical and Alignment Hardware

Optical bench with graduated rails

Lens holders and mounts

XYZ translation stage

Rotary stage for angular characterization

5. Electrical Measurement Instruments

Source-measure unit (SMU) for I-V curve measurement

Digital storage oscilloscope

DC regulated power supply

Picoammeter or electrometer (for low currents)

Digital multimeter ($\mu\text{A}/\text{mA}$ ranges)

Function generator for modulated light studies

6. Optical Power Measurement

Optical power meter with calibrated sensor head

Spectroradiometer or optical spectrum analyzer (for wavelength response)



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Apparatus for Opto Electronics Characterization

KS/PHY/OPTOELECHAR

Characteristic study of Light Dependent Resistor (LDR), Characteristic study of Light Emitting Diode (LED), Characteristic study of Photo Transistor, Characteristic study of Photo Diode, Characteristic study of Solar Cell, Characteristics of Opto - Coupler

Required Components For the experimental Set-Up

Complete Experimental Set up

Required Components For the experimental Set up

Optical Rail

Length: 500mm

Material: Aluminum alloy

Finish: Black anodized

Quantity: 1 no

Detector Mount

Material: Black anodized aluminium alloy

Diameter: 30mm

Quantity: 1 no

Diode Laser with Power Supply (Red)

Wavelength: 650 nm

Output Power: 3 mW

Quantity: 1 no

Kinetic Laser Mount

Material: Black anodized aluminium alloy

Adjustment: Using 80 tpi lead screws

Adjustment Range: +/- 4 degrees

Quantity: 1 no



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Apparatus for Opto Electronics Characterization

KS/PHY/OPTOELECHAR

Mounted Opto-Electronic Detectors

Light dependent resistor (LDR)

Light emitting diode (LED)

Photo Transistor

Photo Diode

Scar Cell

Opto-Coupler

Quantity: 1 no each

Polarizer Rotator with Mount

Material: Sheet Polarizer

Rotation: 360 degree

Quantity: 1 no

Resolution : 1degree Quantity

Opto Electronic measurement Unit

Output Voltage

Port 1: 0-5 Volt Variable

Part 2: 5 Volt fixed

Quantity: 1 no

Pc is not scope of supply.



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com



Carey Foster Bridge Apparatus

KS/PHY/CAREYFBRDG

Purpose

The Carey Foster Bridge is used to measure very small differences in resistance with high precision. It is especially useful for calibrating low-value resistors.

Principle

It is based on the Wheatstone bridge principle. Instead of directly measuring an unknown resistance, the bridge detects the *difference* between two nearly equal resistances by shifting the balance point along a uniform slide wire.

Apparatus Components

Carey Foster bridge base with uniform slide wire (typically 1 m).

Set of known standard resistances.

Unknown low-value resistor.

Two nearly equal ratio arms P and Q

Galvanometer.

Battery or regulated DC supply.

Key/switch and connecting wires.

Determination of unknown resistance by Carey Foster Bridge method

Required Components For the experimental Set-Up

Complete Experimental Set up

Required Components For the experimental Set-Up

Meter bridge, unknown resistance (5 nos), variable rest box (2 nos), DC power supply (2nos), fixed resistance(2 nos), connecting wire, extra bridge wire



Kalpana Scientific Pvt.Ltd.

Konnagar, West Bengal-712235, India

Tel. 7718610048 website: www.kalpanascientific.com

Email id : kalpanascientifickolkata@gmail.com

